



TUBE FEEDING ADMINISTRATION VIA ENTERALITE INFINITY PUMP

GETTING STARTED

- Use a clean work area.
- Clean your hands thoroughly.

GATHER YOUR SUPPLIES

- Formula (room temperature)
- EnteraLite Infinity pump with charging cable
- Administration set (feeding bag with tubing)
- Backpack

PROGRAM THE PUMP

- Press and hold ON/OFF key for 1.5 seconds to turn pump on.
- Program the rate. Press RATE/DOSE key to display rate. Use the (+) and (-) keys to set pump to desired rate.
*Rate is the amount of formula infused per hour measured in mL/hr.
- Program the dose. Press RATE/DOSE key to display dose. Use the (+) and (-) keys to set pump to desired rate.
 - Unless otherwise noted, press (+) key until pump reaches infinity (INF) - INF is reached after exceeding 3000mL. This setting will enable the pump to run any time there is formula in the feeding bag.
- FEED INT key will program interval feeding. This setting does not need to be programmed unless otherwise noted by educator or physician.
- CLEAR key will clear out any setting you had previously entered.
- VOL/TOTAL key will display volume delivered for current feeding. Press key again to display total cumulative volume delivered in all feed cycles since total volume was last cleared.
- Keep the pump plugged in to the power cord when able as battery life will only last 24 hours (if running at dose of 125ml/hr or less). Align ridges on power cord with ridges on pump to ensure accurate connection.

PREPARE FOR FORMULA ADMINISTRATION

- Pour formula into feeding bag. Remove as much air as possible and close the lid by pushing down and turning.
- Load the tubing. Open the pump door and put the tubing loop over the pump wheel. Pull the tubing to the right until the cassette fits into the pump correctly. Close the pump door.
- Prime tubing. Turn on pump by pressing ON/OFF key for 1.5 seconds. Hold down the PRIME key on your pump to fill the tubing with formula and prevent air from entering the GI tract.



- Place pump and administration set into backpack or hang on IV pole. The back pocket is for the feeding bag. The small center pocket is for an ice pack if needed. The front pocket is for the pump and there is a small front pocket for easy viewing of pump screen.

FORMULA ADMINISTRATION

- Attach the end of the tubing to your feeding tube. (If you have an EnFit feeding tube, screw off tip of tubing to expose EnFit adapter).
- Start the pump by pressing the RUN/PAUSE key. If the pump is running, the screen will show the word run with a circle flashing around it.
- If an alarm displays while the pump is running, check that the tubing does not have any air bubbles, kinks, or bends. Turn the pump off and start it again. Check the troubleshooting manual if this does not resolve the issue.

HELPFUL TIPS

- Store unopened formula at room temperature, always check expiration date before administration.
- A maximum 12-hour hang time is recommended for formula, and the administration set should be flushed with water before adding additional formula.
- Administration set should be changed every 24 hours.
- Cover open cans of formula and store them in the refrigerator between feedings. Throw away any open, unused cans of formula after 24 hours.
- To avoid tube clogging, always flush your feeding tube with the amount of water recommended by your healthcare provider (at least 30mL) before and after:
 - Each feeding.
 - Taking each medication.
- If you are not using your feeding tube daily, flush it with 60 mL of water at least once a day to avoid clogging.
- IF your feeding tube gets clogged, flush with 60mL of warm water and clamp tube to let water sit in tube and break down the clog. If this does not unclog the tube, call your healthcare provider.
- Move the clamp to different sections of your feeding tube to avoid permanent kinking the tube.

Irving: 8200 Tristar Drive, Suite 120
Irving, Texas 75063

Fort Worth: 3515 NW Jim Wright Freeway
Fort Worth, Texas 76106

Pharmacy Phone: (866) 388-3883
Pharmacy Fax: (888) 249-3132

www.landmarkinfusion.com